

# grid\_scan

March 15, 2021

## 1 Grid scan and time series at the beamline

### 1.1 Load the name space of the beamtime

```
[1]: # load the beamtime name space
      %run -i 999-load2.py
```

```
/Users/sst/anaconda3/envs/test_xpystack/lib/python3.8/site-
packages/databroker/v1.py:1602: UserWarning: Failed to load config. Falling back
to v0.Exception was: Unable to handle metadatastore.module
'databroker.headersource.sqlite'
      warnings.warn(
```

No config file could be found in the following locations:

```
/Users/sst/.config/acq
/Users/sst/anaconda3/envs/test_xpystack/etc/acq
/etc/acq
```

Loading from packaged simulation configuration

INFO: Initializing the XPD data acquisition environment ...

INFO: area detector has been configured to new acquisition time (time per frame)  
= 0.1s

INFO: Reload beamtime objects:

ScanPlans:

```
0: ct_5
1: ct_0.1
2: ct_1
3: ct_10
4: ct_30
5: ct_60
6: Tramp_5_300_310_3
```

Samples:

```
0: Ni
```

```
{'Verification time': '2021-03-15 20:17:04', 'Verified by': 'st'}
```

Is this configuration correct? y/n: y

Please input your initials: st

INFO: beamtime object has been linked

INFO: Initialized glbl, bt, xrun.

INFO: Publish data to localhost port 5567 with prefix 'raw'.

INFO: Changed home to /Users/sst/acqsim/xpdUser

OK, ready to go. To continue, follow the steps in the xpdAcqdocumentation at <http://xpdacq.github.io/xpdacq>

## 1.2 Import the packages we need

```
[2]: import scanplans.planhelper as pth
import bluesky.plan_stubs as bps
import bluesky.plans as bp
import matplotlib.pyplot as plt
```

Here, I use some fake motors and a fake two dimensional detector.

```
[3]: # fake devices
from ophyd.sim import hw
HW = hw()
MOTORX = HW.motor1
MOTORY = HW.motor2
MOTORPHI = HW.motor3
DETECTOR = xpd_configuration['area_det']
```

## 1.3 Set the exposure time

```
[4]: plan = pth.configure_area_det(DETECTOR, 0.1, 0.1)
xrun({}, plan)
```

INFO: requested exposure time = 0.1 - > computed exposure time= 0.1

```
[4]: ()
```

## 1.4 Grid scan


### 1.4.1 A one dimensional scan

```
[5]: plan5 = pth.xpdacq_grid_scan([DETECTOR], MOTORX, -1, 1, 3)
uids5 = xrun(0, plan5)
```

```
INFO: closing shutter...
INFO: taking dark frame...
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
dark frame complete, update dark dict
opening shutter...
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

### 1.4.2 Visualize the motor position and images in a scan

```
[6]: data = db[uids5[-1]].primary.read()
data["motor1"].plot(x="time", marker="o")
plt.show()
```




### 1.4.3 A one dimensional grid list scan

```
[7]: plan7 = pth.xpdacq_list_grid_scan([DETECTOR], MOTORPHI, [0, 15, 45])
      uids7 = xrun(0, plan7)
```

```
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

### 1.4.4 Visualize the position of the motor

```
[8]: data = db[uids7[-1]].primary.read()
      data["motor3"].plot(x="time", marker="o")
      plt.show()
```



#### 1.4.5 A two dimensional grid scan


We scan the sample by moving the MOTORX and MOTORY. The MOTORX and MOTORY will move from -1 to 1 and 3 points will be measured on the grid so there are 3 \* 3 images in total.

```
[9]: plan0 = pth.xpdacq_grid_scan([DETECTOR], MOTORY, -1, 1, 3, MOTORX, -1, 1, 3)
      uids0 = xrun(0, plan0)
```

```
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

#### 1.4.6 Visualize the motor movement in a grid scan

```
[10]: data = db[uids0[-1]].primary.read()
      data.plot.scatter("motor1", "motor2", hue="time", size=5, aspect=1.2)
      plt.plot(data["motor1"], data["motor2"], '--')
      plt.show()
```




#### 1.4.7 A two dimensional scan where a dark is taken at the point that the slow motor changes its position

```
[11]: plan6 = pth.xpdacq_grid_scan_with_dark([DETECTOR], MOTORY, -1, 1, 3, MOTORX,
↪ -1, 1, 3)
uids6 = xrun({}, plan6)
```

```
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

#### 1.4.8 Visualize the motor movement and where the dark is taken

```
[12]: data = db[uids6[-1]].primary.read()
darkData = db[uids6[-1]].dark.read()
facet = data.plot.scatter("motor1", "motor2", hue="time", size=5, aspect=1.2)
facet.axes.scatter(darkData["motor1"], darkData["motor2"], marker="x", c="r",
    ↪label="where to take dark")
facet.axes.legend()
plt.plot(data["motor1"], data["motor2"], '--')
plt.show()
```



#### 1.4.9 A two dimensional snake scan


```
[13]: plan1 = pth.xpdacq_grid_scan([DETECTOR], MOTORY, -1, 1, 3, MOTORX, -1, 1, 3,
    ↪snake_axes=True)
uids1 = xrun(0, plan1)
```

INFO: This scan will append calibration parameters recorded in  
/Users/sst/acqsim/xpdUser/config\_base/xpdAcq\_calib\_info.poni

```
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

#### 1.4.10 Visualize the motor movement in a snake scan

```
[14]: data = db[uids1[-1]].primary.read()
data.plot.scatter("motor1", "motor2", hue="time", size=5, aspect=1.2)
plt.plot(data["motor1"], data["motor2"], '--')
plt.show()
```






#### 1.4.11 A three dimension scan

```
[4]: plan2 = pth.xpdacq_grid_scan([DETECTOR], MOTORY, -1, 1, 3, MOTORX, -1, 1, 3,
    ↪MOTORPHI, -30, 30, 3, snake_axes=[MOTORX])
uids2 = xrun(0, plan2)
```


```
INFO: closing shutter...
INFO: taking dark frame...
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
dark frame complete, update dark dict
opening shutter...
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

#### 1.4.12 Visualize the motor movement

```
[8]: data = db[uids2[-1]].primary.read()
data.plot.scatter("motor1", "motor2", hue="time", size=5, aspect=1.2)
plt.plot(data["motor1"], data["motor2"], '--')
plt.show()
```



```
[9]: fig = plt.figure(figsize=(8, 8))
ax = plt.axes(projection='3d')
ax.plot3D(data["motor1"], data["motor2"], data["motor3"], '--')
scatter = ax.scatter3D(data["motor1"], data["motor2"], data["motor3"],
    ↪c=data["time"], s=60)
ax.set_xlabel('motor1')
ax.set_ylabel('motor2')
ax.set_zlabel('motor3')
plt.show()
```



## 1.5 Time series

### 1.5.1 Configure the velocity of a motor

```
[18]: plan = bps.mv(MOTORPHI.velocity, 10)
      xrun({}, plan)
```

[18]: ()

### 1.5.2 A helper function to calculate the time and positions in a fly scan

```
[19]: velocity = pth.calc_velocity(start=0, end=30, exposure=10, num=2)
      plan = bps.mv(MOTORPHI.velocity, velocity)
      xrun({}, plan)
```

```
[19]: ()
```


### 1.5.3 A time series of one detector and one motor (use ramp plan)

```
[20]: # move the motor to 30, trigger and read dectector until it finishes
      plan3 = pth.xpdacq_ramp_count([DETECTOR, MOTORPHI], MOTORPHI, 30)
      uids3 = xrun(0, plan3)
```

```
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

### 1.5.4 Visualize the data

```
[21]: data = db[uids3[-1]].primary.read()
      data["motor3"].plot(x="time", marker='o')
      plt.show()
```



### 1.5.5 A time series (use set and count)


```
[22]: # move the motor to 30 and trigger and read dectector for 4 times in total and
      ↪ the period between the start of each trigger is 1 s.
      plan4 = [bps.abs_set(MOTORPHI, 30, wait=False), pth.xpdacq_count([DETECTOR,
      ↪ MOTORPHI], 4, 1)]
      uids4 = xrun(0, plan4)
```

```
INFO: This scan will append calibration parameters recorded in
/Users/sst/acqsim/xpdUser/config_base/xpdAcq_calib_info.poni
INFO: Current filter status
INFO: flt1 : In
INFO: flt2 : Out
INFO: flt3 : Out
INFO: flt4 : Out
```

### 1.5.6 Time series of the motor position

```
[23]: data = db[uids4[-1]].primary.read()
      data["motor3"].plot(x="time", marker='o')
```

```
[23]: [<matplotlib.lines.Line2D at 0x7ffb9a1b0a30>]
```



## 1.6 Visualization of images

Visualize the image in an array using the jupyter widget.

```
[12]: data = db[uids2[-1]].primary.read()
```

```
def visImage(index, **kwargs):
    data["pe1_image"][index].plot.imshow(aspect=1.2, size=5, **kwargs)
    plt.show()
```

```
[22]: from ipywidgets import interact
```

```
interface = interact(visImage, index=(0, 26), vmax=(0., 3e4), vmin=(0., 3e4))
```

```
interactive(children=(IntSlider(value=13, description='index', max=26), FloatSlider(value=15000, description='vmax', max=30000)))
```

```
[ ]:
```